

Summary of the most commonly used test methods

Part 13 – Needle-based Injection Systems (ISO 11608-1)

The requirements and test methods for needle-based injection systems are described in ISO 11608-1.


□ Pre-conditioning




Cold storage




Dry heat



Damp heat



Vibration





Free fall testing

□ Dose-accuracy testing

Last-dose testing

Life-cycle testing



For further information also regarding additional test methods not covered by this summary please contact:

SpineServ GmbH & Co. KG, Soefflinger Strasse 100, 89077 Ulm, Germany

Phone: +49 (0)731 175 6788

mail@spineserv.de

www.spineserv.de