



Summary of the most commonly used test methods

Stainless steel needle tubings


ISO 9626 - Stainless steel needle tubing for the manufacture of medical devices — Requirements and test methods



ISO 9626, Chapter 5.8 (Annex B): Stiffness of tubing



ISO 9626, Chapter 5.9 (Annex C): Resistance of tubing to breakage



ISO 9626, Chapter 5.10 (Annex D): Resistance to corrosion

For further information please contact:

SpineServ GmbH & Co. KG, Soefflinger Strasse 100, 89077 Ulm, Germany

Phone: +49 (0)731 175 6788

experts@spineserv.de

www.spineserv.com